

DAKEWE

SurePrint C7 cassette printer

The compact allrounder



SurePrint C7 cassette printer

The future is non-contact laser marking

Enjoy smooth printing performance and savings on operational expenditure provided by non-contact UV laser printing technology

Compact and well-packed

Fits neatly into grossing station and labs

Continuous output

Special-engineered cassette output rack allows for non-stop printing, never stop the workflow

Omnipotent to print on all sorts of cassettes

Capable to print all regular cassette configurations: Front or back hinged, with or without lid, one-piece or separate

Intuitive handling and viewing of printed cassettes

UV laser tackles a wide range of polymers and prints well on different cassettes



Simple and powerful

Simple operation with friendly UI, yet powerful to adapt to various lab requirements

Easy integration with various LIS

Flexible template editing to provide layout fulfilling every need

Rich content including multiple barcode formats

Fast printing, long-lasting marks

Great barcode scannability thanks to 2,500 dpi high resolution

Permanent marks that resist processing solvents, heat and scratch

Hassle free maintenance

No more ink cartridges, ribbons, light bulbs or print heads

Eliminates the costs and downtime associating with these consumables

Care for the operators

Dedicated air purification system keeps the environment fresh

Smooth touch on all edges and corners



Specifications

Printing technology	Non-contact UV laser marking technology
Print content	Numbers, characters in multiple languages, symbols, 1D/2D barcodes, graphics
Print speed	Up to 2.5 sec/cassette
Cassette type	Cassettes with/without lid, one-piece cassettes; 35° - 45° angle
Air purification	Build-in air purification system with carbon filter
Load capacity	Up to 80 cassettes in 1 hopper
Output rack	Continuous output
System interface	LIS, HIS, PIC, PACS
Operating systems	Windows 7, Windows 10, Windows 11, Linux
Power	100 - 240 VAC, 50/60 Hz, 210 VA
Dimensions	320 × 380 × 440 mm
Weight	Approx. 25 kg

DAKEWE

Dakewe medical
Shenzhen, China
www.dakewemedical.com
Ver 1.0, Rev. May 2024



CasaÁlvarez

MADRID: T. 91 548 10 16
atencionalcliente@casaalvarez.com
BARCELONA: T. 93 433 44 00
cabarna@casaalvarez.com
GRANADA: T. 679 59 69 12
ca.granada@casaalvarez.com
VALENCIA: T. 634 95 20 79
ca.valencia@casaalvarez.com
SEVILLA: T. 683 14 40 56
ca.sevilla@casaalvarez.com

www.casaalvarez.com

